



INSTRUCTION

WorkBeads 40 TREN

WorkBeads™ 40 TREN resin for multimodal ion exchange chromatography (IEX) has a ligand that is positively charged below approx. pH 9. This resin can be used for several different applications, especially due to its higher salt tolerant properties, e.g., for alternative IEX selectivity, for sample cleanup in monoclonal antibody (mAb) purification processes to guard the protein A column from viruses and different host cell impurities, or as a polishing step in the mAb purification process

- Differential selectivity due to higher salt tolerance and multimodal properties
- Reduced fouling of e.g., protein A resins by viruses and different host cell impurity removal
- High binding capacity and purity



Intended use

WorkBeads resins are developed and supported for both research and production-scale chromatography. WorkBeads resins are produced according to ISO 9001:2015, and Regulatory Support Files (RSF) are available to assist the process validation and submissions to regulatory authorities.

GoBio™ prepacked column family is developed for convenient, reproducible, and fast results and can be used from small scale purification through process development to full-scale manufacturing.

Safety

Please read the associated Safety Data Sheet (SDS) for WorkBeads 40 TREN and the safety instructions for any equipment to be used.

Unpacking and inspection

Unpack the shipment as soon as it arrives and inspect it for damage. Promptly report any damage or discrepancies to complaints@bio-works.com

Short protocol

This short protocol is for column packing of WorkBeads 40 TREN resin and general steps for purification. Detailed instructions and recommendations for optimization are given later in this instruction.

1. Make a slurry of the desired resin concentration.
2. Pour the slurry into the column.
3. Pack the resin with an appropriate flow rate.
4. Apply an axial compression of less than 2%.
5. Equilibrate the column with binding buffer.
6. Apply sample.
7. After sample application, remove unbound material by washing with, e.g., 20 – 30 CV (column volumes) washing buffer.

Note: If the resin is used in “negative mode” the target molecule will elute in the wash step.

8. Elute the target protein with elution buffer.
9. Wash the column with deionized water.
10. Equilibrate the column with 20% ethanol for storage.

Principle

WorkBeads 40 TREN resins contain ligands based on Tris(2-aminoethyl)amine (TAEA). The structure of the ligand used in WorkBeads 40 TREN is shown in Figure 1.




Figure 1. Structure of the ligand used in WorkBeads 40 TREN.

Multimodal ion exchange chromatography

Multimodal ion exchange chromatography separates protein, peptides and other biomolecules via a ligand acting with more than one interaction sites. The interaction utilizes two, or more, different properties for example charge and hydrophobicity. Depending on the chromatographic conditions the interactions differ and work either together or separate in the purification procedure.

Using multimodal ion exchange chromatography can be an excellent alternative, when the common techniques (e.g., ion exchange chromatography, size exclusion chromatography or affinity chromatography) are insufficient. However, to reach the optimal purification using multimodal chromatography, the purification process needs to be optimized according to the properties of the target molecule.

WorkBeads 40 TREN in mAb purification

Purification of monoclonal antibodies usually involves purification on Protein A chromatography resins followed by polishing steps based on anion- or cation exchange chromatography.

The use of WorkBeads 40 TREN in binding or flow through mode will also facilitate removal of nucleic acids, viruses, host cell proteins and other cell-derived impurities. As protein A ligands may be cleaved by proteases, leached protein A ligands can be removed by a polishing step using WorkBeads 40 TREN after the protein A purification step. Notice that the majority of mAbs are basic, thus are mainly positively charged at neutral pH or low pH, and therefore do not bind to the resin.

The characteristics of WorkBeads 40 TREN can be exploited in several ways in mAb purification process:

1. As a guard column for removal of different host cell impurities before the protein A column.
2. In a polishing step after the protein A purification step.

Column packing

Columns with i.d. ≤ 10 mm

WorkBeads resins are cross-linked using a proprietary method that results in a very rigid resin that tolerate pressures of several bars and consequently can run at high flow rates. Follow both these general advices when packing a column and the column manufacturer's specific instructions. Preferably, use a column with an adjustable adaptor. In some instances, a packing reservoir or column extension may be needed.

Note: Always make sure that the maximum pressure of the column hardware is not exceeded. Backpressure caused by the chromatography system components connected downstream of the column may reduce the maximum flow that can be used. Wear eye protection.

For laboratory scale columns, we recommend a so-called single step flow packing procedure since the contribution of the wall support is significant. A single step flow packing is fast and easy to implement. The wall support phenomena are described in the next section.

1. Wash the resin

Most WorkBeads resins are supplied in 20% ethanol. To avoid undue backpressure when packing, wash the desired amount of resin with several column volumes (CV) of deionized water or packing buffer before packing.

2. Make a slurry

Add packing solution to the washed resin to obtain a 40% to 60% slurry concentration. Make approximately 15% extra slurry to compensate for the compression during packing. The amount of slurry to be prepared can be calculated according to:

$$\text{Slurry volume} = \frac{\text{bed volume} \times 100}{\% \text{ slurry}} \times 1.15$$

The concentration of slurry can be determined using different methods. One such method is the centrifugation method: Transfer homogenized slurry to a graded centrifuge test tube. Centrifuge for 3 minutes at 3000 rpm and determine the relation between total volume and resin volume. If the relation is not within desired range, modify the slurry concentration of the resin to be packed by either adding or removing packing solution.

Note: If there is a limitation on the total slurry volume, slurry concentrations up to 70% can be used to compensate the lower volume.

Note: Different packing buffers can be used, such as 0.4 M NaCl or 20% ethanol. Some resins are better packed in packing buffer with higher conductivity.

3. Pour the slurry into the column

Pour the slurry slowly down the side of the column to avoid formation of air bubbles. Preferably, use a packing reservoir or an extra column tube to extend the column volume to accommodate the entire slurry volume during packing. If no packing adaptor is available, packing can be done by stepwise additions and packing. Although not recommended this will give acceptable results for most applications.

4. Pack the bed

Pack the resin with a downward flow higher than the intended operational flow. We recommend 600 cm/h for columns up to 10 mm i.d. Make sure the packing flow rate does not exceed the maximum pressure of the column hardware or the resin. The operational flow should not be more than 75% of the packing flow rate.

$$\text{Volumetric flow rate (ml/min)} = \frac{\text{Linear flow (cm/h)}}{60} \times \text{column cross sectional area (cm}^2\text{)}$$

$$\text{Column cross sectional area (cm}^2\text{)} = \frac{\pi \times d^2}{4}$$

d = column inner diameter in cm

After flow compression, and before mechanical compression, the bed height can be adjusted (i.e., using a spatula) so that the exact bed height is achieved, but do this with care since it is important not to disturb the flow compressed bed left in the column.

5. Close the column

When the bed height is constant mark the bed height on the column. Stop the flow and remove the packing reservoir or extra tube. Note that the bed height will increase temporarily when the flow is stopped. If needed, adjust the bed height by removing excess resin but be careful not to remove too much. Gently fill the column with packing solution to its rim without disrupting the packed bed. Insert the adjustable adapter on top of the packed bed. Apply an axial compression of 1 – 2% of the final bed height by lowering the adapter into the packed bed.

6. Apply a flow (conditioning of bed)

Apply a flow of 450 cm/h and check for any gap formation above the surface of the resin bed. If a gap is observed, stop the flow and adjust the adaptor to eliminate the gap.

Columns with i.d. \geq 10 mm

During scale up, the phenomena of wall support gradually is decreased/diminished. Wall support means that the column wall supports the resin bed and gives it better flow properties. This mainly occurs when the column inner diameter is smaller, i.e. \leq 10 mm. To accommodate for the lack of wall support in wider columns, we recommend a two-step flow packing procedure. The bed should first be consolidated at a low flow rate to ensure the bed is optimally settled (60 cm/h for at least 1.5 CV). After the initial consolidation step, the bed is further compressed at a higher flow rate (300 – 450 cm/h depending on restrictions for at least 1.5 CV). When the bed has been flow-packed, a mechanical compression of 10 – 20% can be performed.

The optimal mechanical compression to apply will depend on the dimensions and resin type, therefore it is essential to evaluate the packed column prior to use (asymmetry and plate number tests as described below).

Note: After flow compression, and before mechanical compression, the bed height can be adjusted (i.e., using a spatula) so that the exact bed height is achieved but this is not to recommend at larger scales.

There is a guide available for large-scale packing,
<https://www.bio-works.com/packing-guide-download>

Evaluation of the packed column

Test the column efficiency to check the quality of packing. Testing should be done after packing, at regular intervals during the usage of the column or when separation performance is seen to deteriorate.

The best method of expressing the efficiency of a packed column is in terms of the height equivalent to a theoretical plate (HETP) and the asymmetry factor (A_s). These values are easily determined by applying a sample such as 1% acetone solution to the column.

For optimal results, the sample volume should be 2.5% of the column volume (CV) and the flow rate 30 cm/h.

If an acceptance limit is defined in relation to column performance, the column plate number can be used as one of the acceptance criteria for column use.

Note: The calculated number of plates will vary according to the test conditions and should only be used as a reference value. Keep test conditions and equipment constant so that results are comparable. Changes of for example solute, solvent, eluent, sample volume, flow velocity, liquid pathway, temperature, etc. will influence the results.

Measuring HETP and A_s

Calculate HETP and A_s from the UV curve (or conductivity curve).

$$\text{HETP} = \frac{L}{N}$$

$$N = 5.54 \times \left(\frac{V_R}{W_h} \right)^2$$

L = bed height (cm)

N = number of theoretical plates

V_R = volume eluted from the start of sample application to the peak maximum

W_h = peak width measured as the width of the recorded peak at half of the peak height

V_R and W_h are in the same units

The concept of reduced plate height is often used for comparing column performance. The reduced plate height, h , is calculated:

$$h = \frac{\text{HETP}}{d_{50v}}$$

d_{50v} = Median particle size of the cumulative volume distribution (cm)

As a guideline, a value of < 3 is very good.

The peak must be symmetrical, and the asymmetry factor as close to 1 as possible. (A typical acceptable range is $0.7 < A_s < 1.3$). A change in the shape of the peak is usually the first indication of bed deterioration.

Peak asymmetry factor calculation:

$$A_s = \frac{a}{b}$$

a = ascending part of the peak width at 10% of peak height

b = descending part of the peak width at 10% of peak height

Figure 2 below shows a UV trace for acetone in a typical test chromatogram from which the HETP and A_s values are calculated.




Figure 2. A typical test chromatogram showing the parameters used for HETP and A_s calculations.

Purification

Multimodal IEX

WorkBeads 40 TREN resin is positively charged below approx. pH 9. Choose a suitable pH and buffer for the binding of the target protein or the impurities (if the purification is done in “negative mode”). One pH unit above pI, for WorkBeads 40 TREN is a good starting point. The binding conditions should be optimized to achieve binding of the target protein or the impurities. When scouting for the best conditions it is important to start with sufficiently low ionic strength. Guideline for a starting point for designing the experiment is given in Table 1.

Table 1. Typical buffer composition for purification using WorkBeads 40 TREN.

Buffer	Buffer composition
Binding buffer	50 mM Tris-HCl, pH 7.4
Elution buffer	50 mM Tris-HCl, 1 M NaCl, pH 7.4

When the target molecule is bound to the resin elution can be carried out by applying a linear gradient of increasing concentration of NaCl, by gradually increasing the proportion of elution buffer (high salt). A short step gradient to 1 or 3 M NaCl for 5 column volumes (CV) can be included after elution to ensure desorption of all interacting proteins. When suitable elution conditions are known it is common to apply step gradient elution.

A cleaning-in-place (CIP) step using 5 CV 1 M NaOH is recommended between all runs and should be followed by a careful re-equilibration before the next run. Make sure that elution with neutral salt, e.g., NaCl, has been done before CIP to avoid precipitation of adsorbed proteins.

Sample preparation

After cell disruption or extraction, clarify the sample by centrifugation at 10 000 – 20 000 × g for 15 – 30 minutes. It is generally also recommended to pass the sample through a 0.22 – 0.45 μm filter, e.g., a syringe filter, to avoid transferring any remaining contaminating particles onto the column. Large sample volumes may be clarified by filtration through depth filters or by tangential flow filtration, which may be cheaper and more efficient than investing in a large-scale centrifuge. Application of a sample that has not been properly clarified may reduce the performance and lifetime of the column.

The sample should have a pH that confers a net charge to the target protein that is opposite to the charge of the column resin if the focus is to bind the target protein. The ionic strength should be low. The optimal binding conditions depends on the combination of the pH and the ionic strength. The sample solution may therefore need to be adjusted before applied to the column. It is generally recommended that the sample should have a similar pH and conductivity as the binding buffer. Sample adjustments can be done by dilution using the binding buffer, by chromatographic desalting or diafiltration, or through adjusting the pH by addition of an acid or base.

When the resin is used as a guard column for removal of chromatin and other impurities before the Protein A column in a mAb purification process the target protein is washed out with the wash buffer after sample application.

Optimization

The following paragraphs will give indications on some parameters that can be tuned to get the optimal conditions for purification of using WorkBeads 40 TREN resin.

Selection of buffer

Selecting a buffer with optimal conditions for the target protein will improve the result of the purification. The buffer should be chosen with a pK_a -value within 0.5 units from the intended pH to obtain a high enough buffering capacity. Table 1 shows one example of buffers which can be used for ion exchange chromatography, however the buffer choice will be depending on the target molecule and aim of the purification procedure. For other useful buffers and their pK_a -values at 25 °C see reference: Methods in Enzymology, Volume 463, pp 46-47, Burgess, R.R and Deutcher M.P.

The buffer substance should be selected to have the same charge as the resin. A buffer with opposite charge will interact with the charged groups in the resin and may cause local pH disturbances that destroys the separation. Usually, low conductivity in the binding buffer is preferred but optimization regards to pH and conductivity can improve binding capacity. An increase in ionic strength may decrease the ability of contaminants to bind while the target protein remains bound. However, chromatographic conditions should be chosen so that the protein is stable during purification.

Optimization of binding conditions

The key conditions to be optimized is usually pH and conductivity (by addition of NaCl or other salts, or dilution). The conditions must also be selected to keep the protein in its native state.

The flow rate during sample loading affects the binding capacity and resolution during the elution. A low flow rate during sample application promotes binding capacity since more time is allowed for mass transport of the target substance into the pores of the resin. A small substance, e.g., a peptide, that has a high diffusion rate will have fast mass transport into the resin and can thus be adsorbed efficiently at high flow rates. A large target substance (e.g., a large protein) has a lower diffusion rate and is more hindered by the walls in the pores giving slow mass transport, causing reduced dynamic binding capacity. A high binding capacity of this substance may thus require a lowered flow rate. If only a part of the binding capacity of the column is used the sample application can be at a higher flow rate without loss of the target substance.

The residence time can be defined as the time between entering and exiting the column of specific part of the sample or buffer. The residence time depend on the flow rate and the dimensions of the column and is typically 1 to 5 minutes in IEX. Typical linear flow rates are 150 – 300 cm/h. See further discussion about flow in the section "Scale-up".

Optimization of washing

A continuously decreasing UV signal is an indication of unbound material still being washed out. The washing should continue until the UV signal is stable and the same as in the washing buffer, or at least not more than 20 mAU. The washing buffer can be the same as the binding buffer, but it may be useful to add an additional step with a dedicated washing buffer to improve purification.

Optimization of elution conditions

Elution can be carried out using a high salt concentration or by altering the pH to change the charge of the adsorbed protein. A stronger binding may require higher salt concentration for elution. The optimal salt concentration is dependent on the purity and recovery requirements as well as the properties of the target protein and the sample. Gradient elution often gives better purity than step elution, but step elution may be desired to obtain the highest possible concentration of the target protein, and to reduce the purification process time. To optimize the salt concentration for step elution an initial gradient test run can be carried out to obtain suitable step elution conditions for purification of the sample, see Figure 3.




Figure 3. Optimization of step elution with salt. A test run with linear gradient elution gives information about suitable salt concentrations to be used in step elution. Note: Remember to take the system dead volume into account when comparing the print outs of the gradient and the trace.

Scale-up

After developing a chromatographic procedure in small scale column, e.g., 7 (i.d.) × 100 mm (bed height), WorkBeads resins can be packed into larger columns for scale-up. Large scale purification is often carried out in columns with bed heights of 200 – 300 mm.

Scale-up principles

During scale-up the ratio between sample volume and column volume should be kept constant. The column volume is scaled up by increasing the column diameter while keeping the bed height the same (e.g., 200 mm). The linear flow rate should remain the same while the volumetric flow rate increases. The volumetric flow rate for each column can be calculated according to:

$$\text{Volumetric flow rate (mL/min)} = \frac{\text{Linear flow rate (cm/h)} \times \text{Column cross sectional area (cm}^2\text{)}}{60}$$

Flow

The concepts of volumetric flow, linear flow rate and residence time are important when scaling-up in chromatography. Volumetric flow is measured in mL/min or L/min, linear flow in cm/h and residence time in minutes. The relationship between these metrics is:

$$\text{Linear flow rate (cm/h)} = \frac{\text{Volumetric flow (mL/min)} \times 60}{\text{Column cross sectional area (cm}^2\text{)}}$$

$$\text{Residence time (minutes)} = \frac{\text{Column bed height (cm)} \times 60}{\text{Linear flow rate (cm/h)}}$$

In the initial process development work it is common to use a small column, e.g., 7 × 100 mm, to save sample, buffers and time. This column has a shorter bed height than the final column which may have a bed height of 200 mm or more. The flow rate for the larger column can be calculated from the flow that was established on the small column, using the equation above by keeping the residence time of the small column the same as for the larger column. This will allow an increase of the linear flow in proportion to the increase in bed height between the columns, see Table 3 for examples. If the column bed height is kept constant during scale-up, the linear flow rate should be kept constant (as well as the residence time).

Table 2. Example of scale-up parameters.

Column dimension	Residence time (minutes)	Linear flow rate (cm/h)	Volumetric flow rate (mL/min)
16x100	4	150	5.0
26x100	4	150	13.3
80x200	8	150	126
130x200	8	150	332
200x200	8	150	785
240x200	8	150	1131
330x250	10	150	2138

Additional purification steps

Optimization of the purification process by tuning the binding, washing and/or elution conditions of the IEX purification step may not be enough to obtain the required purity. Combining two or more purification step based on additional chromatography techniques is then recommended. For example, cation exchange chromatography and anion exchange chromatography can be combined in a purification process. Other techniques, such as size exclusion chromatography (gel filtration) and hydrophobic interaction chromatography (HIC) are commonly used alternatives. Each purification step should be thoroughly optimized, and preferably in the context of the other steps applied on the overall process.

Desalting and buffer exchange

Buffer exchange or desalting of a sample can be used before analysis and/or after purification. This can be carried out quickly and easily in lab-scale using GoBio Mini Dsalt 1 mL, GoBio Mini Dsalt 5 mL, GoBio Prep 16x100 Dsalt (20 mL) and GoBio Prep 26x100 Dsalt (53 mL) prepacked columns depending on sample volumes. GoBio Prod prepacked columns starting from 1 L are available for larger sample volumes, see "Related products".

These columns are very useful alternatives to dialysis for larger sample volumes or when samples need to be processed rapidly to avoid degradation.

Pre-swollen WorkBeads Dsalt is also available in bulk for packing column format of choice.

To find out more about Bio-Works' chromatography products visit www.bio-works.com

Maintenance of the resin

Cleaning and sanitization

During purification impurities such as cell debris, lipids, nucleic acids and protein precipitates from the samples may gradually build-up of contaminants in the resin, or fouling. The severity of this process depends on the type of sample applied to the column, and the pre-treatment of the sample. The impurities covering the resin may reduce the performance of the column over time.

Regular cleaning (cleaning-in-place, CIP) keeps the resin clean, reduces the rate of further fouling, and prolongs the capacity, resolution and flow properties of the column. Cleaning of a column using 1 M NaOH applied by a low reversed flow for 2 hours or overnight is often sufficient.

According to internal studies WorkBeads 40 TREN has very good tolerance towards commonly used CIP solutions (168 h at 40°C).

Sanitization (reduction of microorganisms) can be carried out using combinations of NaOH and ethanol, e.g., incubation with a mixture of 0.5 M NaOH and 40% ethanol for 3 hours. The sanitization procedure and its effectiveness will depend on the microorganisms to be removed and needs to be evaluated for each case.

Storage

Store the resin at 2 to 25°C in 20% ethanol.

Note: Use a reduced flow rate during equilibration with 20% ethanol, maximum 50% of the maximum flow rate.

Product information

WorkBeads 40 TREN	
Target substances	Proteins, peptides, oligonucleotides, viruses
Matrix	Rigid, highly cross-linked agarose
Average particle size (D_{v50}) ¹	45 µm
Ligand	Tris(2-ethylaminoethyl)amine (TAEA)
Ionic capacity	130 – 200 µmol Cl ⁻ /mL resin
Dynamic binding capacity	50 mg BSA/mL resin ²
Max flow rate (20 cm bed height and 5 bar)	600 cm/h
Chemical stability	Compatible with all standard aqueous buffers used for protein purification. Should not be stored at low pH for prolonged time.
Operational pH range ³	2 – 13
CIP and screening pH range ³	2 – 14
Storage	2 to 25 °C in 20% ethanol

¹ The median particle size of the cumulative volume distribution.

² Dynamic binding capacity determined at 4 minutes residence time (0.25 mL/min in 1 mL column) in 50 mM Tris-HCl, 50 mM NaCl, pH 8.0. Optimal flow rate during binding is depending on the sample.

³ Within the operational pH range, the resin can be operated without significant change in function. Within the CIP (Cleaning-in-place) and screening pH range the resin can be subjected to the denoted pH range without significant change in function.

GoBio prepacked column family

GoBio prepacked column family is developed for convenient, reproducible and fast results and includes columns with different sizes and formats.

GoBio Mini 1 mL and GoBio Mini 5 mL for small scale purification and screening using a shorter packed bed.

GoBio Screen 7x100 (3.8 mL) for reproducible process development including fast and easy optimization of methods and parameters.

GoBio Prep 16x100 (20 mL) and GoBio Prep 26x100 (53 mL) for lab-scale purifications and scaling up.

GoBio Prep 16x600 (120 mL) and GoBio Prep 26x600 (320 mL) for preparative lab-scale size exclusion chromatography.

GoBio Prod 80x200 (1 L), GoBio Prod 130x200 (2.7 L), GoBio Prod 200x200 (6 L), GoBio Prod 240x200 (9 L) and GoBio Prod 330x250 (21.4 L) for production-scale purifications.

Related products

Product name	Pack size ¹	Article number
Prepacked columns		
GoBio Mini TREN 1 mL	1 mL × 1	45 655 211
	1 mL × 5	45 655 213
	1 mL × 10	45 655 214
GoBio Mini TREN 5 mL	5 mL × 1	45 655 215
	5 mL × 5	45 655 217
	5 mL × 10	45 655 218
GoBio Mini IEX Screening kit ²	1 mL × 4	45 900 001
GoBio Mini Dsalt 5 mL	5 mL × 5	45 360 107
GoBio Screen 7x100 40 TREN	3.8 mL × 1	55 4630 01
GoBio Prep 16x100 40 TREN	20 mL × 1	55 463 021
GoBio Prep 16x100 Dsalt ³	20 mL × 1	55 700 021
GoBio Prep 26x100 40 TREN ³	53 mL × 1	55 463 031
GoBio Prep 26x100 Dsalt	53 mL × 1	55 700 031
GoBio Prod 80x200 40 TREN ³	1 L × 1	55 463 042
GoBio Prep 80x200 Dsalt ³	1 L × 1	55 700 042
GoBio Prod 130x200 40 TREN ³	2.7 L × 1	55 463 062
GoBio Prep 130x200 Dsalt ³	2.7 L × 1	55 700 062
GoBio Prod 200x200 40 TREN ³	6 L × 1	55 463 072
GoBio Prep 200x200 Dsalt ³	6 L × 1	55 700 072
GoBio Prod 240x200 40 TREN ³	9 L × 1	55 463 082
GoBio Prep 240x200 Dsalt ³	9 L × 1	55 700 082
GoBio Prod 330x250 40 TREN ³	21.4 L × 1	55 463 093
GoBio Prep 330x250 Dsalt ³	21.4 L × 1	55 700 093
Bulk resins		
WorkBeads Dsalt	300 mL	40 360 003
	1 L	40 360 010
	5 L	40 360 050
	10 L	40 360 060

¹ All different pack sizes are available on www.bio-works.com

² GoBio Mini IEX Screening Kit includes one of each: GoBio Mini S 1 mL, GoBio Mini Q 1 mL, GoBio Mini DEAE 1 mL and GoBio Mini TREN 1 mL.

³ Packed on request.

Ordering information

Product name	Pack size	Article number
WorkBeads 40 TREN	25 mL	40 603 001
	150 mL	40 603 003
	1 L	40 603 010
	5 L	40 603 050
	10 L	40 603 060

Orders: sales@bio-works.com or contact your local distributor.

For more information about local distributor and products visit www.bio-works.com or contact us at info@bio-works.com

bio-works.com

Bio-Works, WorkBeads and GoBio are trademarks of Bio-Works Technologies.
All third-party trademarks are the property of their respective owners.

© Bio-Works.

All goods and services are sold subject to Bio-Works terms and conditions of sale.
Contact your local Bio-Works representative for the most current information.

Bio-Works, Virdings allé 18, 754 50 Uppsala, Sweden. For local office contact information, visit bio-works.com/contact.

IN 40 600 020 BB