


# Optimization and scale-up of oligonucleotide purifications



Cecilia Unoson<sup>1</sup>, Lars Holmberg<sup>2</sup>, Kajsa Eriksson<sup>1</sup>, Anna Heijbel<sup>1</sup> and Lars Haneskog<sup>1</sup>

<sup>1</sup> Bio-Works, Uppsala, Sweden, <sup>2</sup> PeptiSystems, Uppsala, Sweden

## Introduction


Solid-phase synthesis of oligonucleotides generally give material of rather high purity. However, for many applications, and especially for therapeutics, there is a need for purification to remove incomplete or erroneous sequences.

Anion exchange chromatography (AEX) is an efficient technique for oligonucleotide purification giving high purity and good yield in a single step, but optimization of the binding and elution conditions is needed. When the process conditions giving required purity and yield have been developed, scale-up can start. We have investigated the change in the purity vs yield correlation when transitioning from low sample load to high sample load.

## Purification optimization

The resolution of oligonucleotide purifications on WorkBeads™ 40Q was compared for NaCl-elution in the presence of 20 mM Tris-HCl, pH 8 and 20 mM NaOH, pH 12 (Fig. 1). The main peak was collected in 1-mL fractions. Each individual fraction was analyzed for yield and purity on a DNAPac PA200 analytical IEX column (Thermo Fisher Scientific). The WorkBeads 40Q resin was also compared to Capto™ Q ImpRes resin (GE Healthcare).

The dynamic binding capacity (DBC) was determined by frontal analysis at 150 cm/h to 48 mg/mL resin by applying the crude oligonucleotide preparation from the solid support without further adjustment of the feed.




**Figure 1.** Purification of oligonucleotides on WorkBeads 40Q using 20 mM NaOH, pH 12 (solid blue) or 20 mM Tris-HCl, pH 8 (dotted green). The elution gradients are shown in red (dotted line 20-40% for 20 mM Tris-HCl, pH 8, and solid line 20-50% for 20 mM NaOH, pH 12).

## Purity and yield

For optimization, selected 1-mL fractions were combined to assess the effect of pooling on yield and purity, see Fig. 2 and Table 1. No significant difference was seen for the two buffer systems used for this oligonucleotide batch. A purity of 95.6% with 74.2% yield was obtained with the broadest pooling using Tris-buffer (Fig. 2B). Higher purities could be obtained by a more narrow pooling, as illustrated in Figure 2B. WorkBeads 40Q gave both higher purity and increased yield compared to Capto Q ImpRes (Table 1).

**Table 1.** Yield and purity. Pools of seven 1 mL fractions of the main peak were compared.

| Resin             | Buffer               | Purity (%) | Yield (%) |
|-------------------|----------------------|------------|-----------|
| Crude preparation | –                    | 84.5       | 100       |
| WorkBeads 40Q     | 20 mM NaOH, pH 12    | 95.5       | 76.9      |
| WorkBeads 40Q     | 20 mM Tris-HCl, pH 8 | 95.6       | 74.2      |
| Capto Q ImpRes    | 20 mM NaOH, pH 12    | 93.7       | 68.9      |
| Capto Q ImpRes    | 20 mM Tris-HCl, pH 8 | 94.4       | 64.4      |




**Figure 2.** (A) Purity (light blue bars) vs yield (green bars) analyzed over the main peak. (B) Purity vs yield plot for different fraction pools (1 mL, 3 mL, 5 mL and 7 mL) of the eluted main peak obtained with 20 mM Tris-HCl, pH 8.

## Scale-up

When the process conditions giving required purity and yield have been obtained, scale-up can be done. To investigate scale-up conditions, sample load of 80 % of the resins DBC, *i.e.*, 132 mg of the oligonucleotide preparation was loaded to the column.

Figure 3A shows a typical recovery shape of the full-length oligonucleotide for a process purification run with a sample load corresponding to 80% of DBC. In the beginning of the elution gradient the N-x species are eluted, and the full-length oligonucleotide starts to elute later in the gradient. The histogram in the chromatogram in Figure 3B shows the purity in individual fractions. This demonstrates the separation of the full-length oligonucleotide.



**Figure 3.** (A) Eluted oligonucleotide peak after loading an amount of crude oligonucleotide corresponding to 80% of the DBC. Green line represent full-length yield measured in each fraction and dotted black line visualizes presence of n-1 and n-x species. (B) Purity of full-length oligonucleotide (light blue bars) was measured in individual fractions. Pools of collected fractions are visualized (as horizontal bars) in the chromatogram.

**Table 2.** Effect on yield and purity using different pooling after the 80% DBC run.

| Pool           | Purity (%) | N-1 (%) | Yield (%) |
|----------------|------------|---------|-----------|
| Pool 1 (55 mL) | 94.5       | 2.1     | 82.2      |
| Pool 2 (45 mL) | 95.2       | 1.8     | 69.9      |
| Pool 3 (35 mL) | 96.3       | 1.6     | 54.3      |

## Conclusions

Purification of a 20-mer oligonucleotide using WorkBeads 40Q gave excellent purity with good yield in both NaOH- and Tris-based buffers. Both purity and yield were higher on WorkBeads 40Q compared to Capto Q ImpRes. The rigidity of WorkBeads 40Q (45-µm bead resin), allows for efficient purifications of full-length oligonucleotides also at process scales.

Capto is a trademark of General Electric Company.