

BabyBio NTA His-tag Screening kit

BabyBio IDA His-tag Screening kit

BabyBio™ His-tag Screening kits contain columns prepacked with WorkBeads™ IDA and WorkBeads NTA charged with Ni²⁺, Co²⁺, Cu²⁺ or Zn²⁺ ions. The kits are excellent tools for screening combinations of metal ions and chelating ligand (NTA or IDA) to optimize purity and yield when purifying polyhistidine-tagged (His-tagged) proteins. Other native proteins containing histidine, cysteine and tryptophan residues may also bind and can therefore be purified using these columns. The selected column can be used to purify up to 70 mg and 350 mg protein respectively using a 1 ml or 5 ml column.

- Pre-charged columns with different metal ions for easy screening for optimal purity
- Ready-to-use columns for fast results
- High binding capacity and purity



Resin description

BabyBio columns provided in the screening kits are prepacked with resins based on cross-linked agarose. The columns are excellent for research scale purification and selectivity screening in process development.

The resins packed in BabyBio His-tag Screening kit columns are WorkBeads 40 NTA or WorkBeads 40 IDA which contain immobilized chelating ligands based on nitrilotriacetic acid (NTA) or iminodiacetic acid (IDA) respectively. The resins are pre-charged with either Ni²⁺, Co²⁺, Cu²⁺ or Zn²⁺ metal ions.

The structures of the chelating ligands used in WorkBeads 40 NTA and WorkBeads 40 IDA are shown in Figure 1.




Figure 1. Structure of the chelating ligands used in WorkBeads 40 NTA (A) and WorkBeads 40 IDA (B) resins.

The main characteristics of BabyBio His-tag Screening kit columns are shown in Table 1. For more details, please see instructions IN 45 700 010.

Table 1. Main characteristics of BabyBio His-tag Screening kit 1 ml and 5 ml columns.

| | BabyBio: Ni-NTA, Co-NTA, Cu-NTA, Zn-NTA | BabyBio: Ni-IDA, Co-IDA, Cu-IDA, Zn-IDA |
|--|--|--|
| Target substance | His-tagged proteins, proteins containing histidine, cysteine and/or tryptophan amino acid side chains | His-tagged proteins, proteins containing histidine, cysteine and/or tryptophan amino acid side chains |
| Resin | WorkBeads 40 Ni-NTA WorkBeads 40 Co-NTA WorkBeads 40 Cu-NTA WorkBeads 40 Zn-NTA | WorkBeads 40 Ni-IDA WorkBeads 40 Co-IDA WorkBeads 40 Cu-IDA WorkBeads 40 Zn-IDA |
| Matrix | Rigid, highly cross-linked agarose | Rigid, highly cross-linked agarose |
| Average particle size ¹ (D _{V50}) | 45 µm | 45 µm |
| Ligand | Nitrilotriacetic acid (NTA) | Iminodiacetic acid (IDA) |
| Metal ion | Ni ²⁺ , Co ²⁺ , Cu ²⁺ or Zn ²⁺ | Ni ²⁺ , Co ²⁺ , Cu ²⁺ or Zn ²⁺ |
| Static binding capacity ² | 70 mg His-tagged protein/ml resin | NA |
| Dynamic binding capacity ² | 50 mg His-tagged protein/ml resin | NA |
| Column volume | 1 ml 5 ml | 1 ml 5 ml |
| Column dimension | 7 x 28 mm (1 ml) 13 x 38 mm (5 ml) | 7 x 28 mm (1 ml) 13 x 38 mm (5 ml) |
| Recommended flow rate | | |
| BabyBio 1 ml | 1 ml/min (150 cm/h) | 1 ml/min (150 cm/h) |
| BabyBio 5 ml | 5 ml/min (225 cm/h) | 5 ml/min (225 cm/h) |
| Maximum flow rate ³ | | |
| BabyBio 1 ml | 5 ml/min (780 cm/h) | 5 ml/min (780 cm/h) |
| BabyBio 5 ml | 20 ml/min (900 cm/h) | 20 ml/min (900 cm/h) |
| Maximum back pressure | 0.3 MPa, 3 bar, 43 psi | 0.3 MPa, 3 bar, 43 psi |
| Chemical stability | Chelating substances (e.g. EDTA) will strip off the metal ions. Stripped column: 10 mM HCl (pH 2), 10 mM NaOH (pH 12), 100 mM sodium citrate-HCl (pH 3), 6 M guanidine-HCl | |
| pH stability | 7 - 9 (working range) 2 - 12 cleaning (stripped resin) Do not keep the resin at low pH for prolonged time. | 7 - 9 (working range) 2 - 12 cleaning (stripped resin) Do not keep the resin at low pH for prolonged time. |
| Storage | 2 to 25°C in 20% ethanol | 2 to 25°C in 20% ethanol |

1. The median particle size of the cumulative volume distribution.

2. The binding capacity is determined using a BabyBio Ni-NTA 1 ml. The binding capacity is dependent on the size of the target protein, and on the competition with other substances.

3. Aqueous buffers at 20°C. Decrease the maximum flow rate if the liquid has a higher viscosity. Higher viscosities can be caused by low temperature (use half of the maximum flow rate at 4°C), or by additives (e.g. use half of the maximum flow rate for 20% ethanol).

Column description

The column is made from biocompatible polypropylene which does not significantly interact with biomolecules. The top and bottom filters are made from polyethylene. The ready-to-use BabyBio columns are delivered with a plug in the inlet, a cut-off outlet and a cap for storage. The columns can be connected to a syringe, pump or chromatography system using finger tight fittings (coned 10–32) for 1/16" o.d. tubing (standard HPLC PEEK tubing).

Applications

BabyBio IMAC columns supplied in the screening kits can easily be used for fast purifications of His-tagged proteins or native proteins containing histidine, cysteine or tryptophan residues. BabyBio Ni-NTA 1 ml and 5 ml columns can be used to purify up to 70 mg or 350 mg of proteins, respectively. Similar capacities can be expected for the other BabyBio IMAC columns. The purity obtained depends on several factors. A sample including impurities that can bind to the resin may reduce the purity of the target protein. Proteins expressed in *E. coli* are usually easier to purify than proteins expressed in eukaryotic systems (e.g., yeast or mammalian cells). The purification result also depends on the structure of the chelating ligand and the nature of the immobilized metal ion. The large selection of BabyBio IMAC columns offers many possibilities, choosing between two different ligands, NTA or IDA, charged with either Ni²⁺, Co²⁺, Cu²⁺ or Zn²⁺ metal ions. BabyBio Ni-NTA is recommended as the starting point for His-tagged protein purification as it in many cases will give excellent purification results. For more difficult purifications, a screening is often recommended with the different BabyBio IMAC columns to find the optimal combination of ligand and metal ion. For very high purity requirements, it is common to add a second purification step to remove the final impurities and for buffer exchange and salt removal. This can be done by using size exclusion chromatography (SEC/gel filtration), such as WorkBeads SEC resins.

Principle

IMAC utilizes the affinity of histidine, cysteine and tryptophan amino acid side chains on the protein surface for binding to transition metal ions, such as Ni²⁺, Co²⁺, Cu²⁺ and Zn²⁺, immobilized via a metal chelating ligand on the chromatography resin.

IMAC is commonly used for purification of recombinant His-tagged proteins. The His-tag is usually composed of six to ten histidyl groups, and is typically placed at the N- or C-terminus of the target protein, although other positions are possible. His-tagged proteins will bind to the chelating ligand (through the metal ion) and the unbound material will pass through the column. The bound proteins are desorbed by stepwise or gradient elution, using a competing agent or lower pH.

Imidazole is recommended for elution. This is the most commonly used competing agent, but histidine, ammonium chloride or histamine can also be used. Before sample application the column should be equilibrated with a low concentration of the competing agent in order to prevent non-specific binding of endogenous proteins that may bind via for example histidine clusters. This is easily done by using the recommended binding buffer.

For more detailed description of the IMAC principle, see instructions IN 45 700 010.

Purification of His-tagged proteins

An example of purification of a recombinant His₆-tagged Green Fluorescent Protein (His₆-GFP) expressed in *E. coli* on eight different pre-charged BabyBio IMAC columns is shown in Figure 2. The purity was checked by SDS-PAGE, shown in Figure 3.

Sample: 2 ml His₆-GFP in binding buffer
 Columns: Pre-charged BabyBio IMAC columns 1 ml
 Binding buffer: 50 mM sodium phosphate, 300 mM NaCl, 10 mM imidazole, pH 8.0
 Elution buffer: 50 mM sodium phosphate, 300 mM NaCl, 300 mM imidazole, pH 8.0
 Gradient: 0 to 100% elution buffer in 15 column volumes (CV)
 Flow rate: 0.5 ml/min (78 cm/h)




Figure 2. Chromatogram showing comparisons of purifications of clarified His₆-GFP on BabyBio NTA 1 ml and BabyBio IDA 1 ml charged with Ni²⁺, Co²⁺, Cu²⁺ and Zn²⁺ ions. The blue and red lines correspond to the absorbance signal at 280 nm and 490 nm, respectively, and the green line to the percentage of elution buffer.




Figure 3. SDS-PAGE analysis of purified His₆-GFP from the previous chromatograms. FT = flow through and EL = eluted protein.

Scale-up

Scale-up can conveniently be carried out from a 1 ml column to a 5 ml column. If increased capacity is required several columns can be coupled in series (column stacking). Note that the backpressure will increase proportionally to the resin bed height (up to a maximum of 5 columns).

Further scale-up can be done by packing bulk WorkBeads IMAC resins in larger columns (see *Related products*). For more detailed description, please see instructions IN 45 700 010.

Cleaning-in-place

During purification impurities such as cell debris, lipids, nucleic acids and protein precipitates from the samples may gradually build up in the resin. The severity of this process depends on the type of sample applied to the column, and the pre-treatment of the sample. The bound impurities may reduce the performance of the packed column over time. Regular cleaning (Cleaning-in-place, CIP) keeps the resin clean, reduces the rate of further contamination, and prolongs the capacity, resolution and flow properties of the column. Cleaning using 1 M NaOH applied by a low reversed flow for 2 hours or overnight is often sufficient.

Before cleaning of IMAC resins the metal ions must be removed from the resin using, for example, 50 mM Na₂EDTA, pH 8.5. After the cleaning, the resin can be re-charged with fresh metal ions.

Sanitization (reduction of microorganisms) can be done using combinations of NaOH and ethanol (e.g., incubation with a mixture of 0.5 M NaOH and 40% ethanol for 3 hours). The sanitization procedure and its effectiveness will depend on the microorganisms to be removed, and needs to be evaluated for each case.

Equipment

Prepacked BabyBio His-tag Screening kit columns can be used with most standard liquid chromatography equipment. Purification can also be carried out using a syringe connected to the column by a luer or a std HPLC connector.

Storage

Equilibrate the columns in 20% ethanol and close it securely using the included plug and cap. Store the column at 2 to 25°C.

Related products

| Product Name | Pack size ¹ | Article number |
|--------------------------|------------------------|----------------|
| Prepacked columns | | |
| BabyBio Dsalt 5 ml | 5 ml x 5 | 45 360 107 |
| BabyBio Ni-NTA 1 ml | 1 ml x 5 | 45 655 103 |
| BabyBio Co-NTA 1 ml | 1 ml x 5 | 45 655 133 |
| BabyBio Cu-NTA 1 ml | 1 ml x 5 | 45 655 123 |
| BabyBio Zn-NTA 1 ml | 1 ml x 5 | 45 655 143 |
| BabyBio Ni-IDA 1 ml | 1 ml x 5 | 45 655 003 |
| BabyBio Co-IDA 1 ml | 1 ml x 5 | 45 655 033 |
| BabyBio Cu-IDA 1 ml | 1 ml x 5 | 45 655 023 |
| BabyBio Zn-IDA 1 ml | 1 ml x 5 | 45 655 043 |
| BabyBio S 5 ml | 5 ml x 5 | 45 200 107 |
| BabyBio Q 5 ml | 5 ml x 5 | 45 100 107 |
| BabyBio DEAE 5 ml | 5 ml x 5 | 45 150 107 |
| OptioBio 40S 10x100 | 7.9 ml x 1 | 55 420 011 |
| OptioBio 40Q 10x100 | 7.9 ml x 1 | 55 410 011 |
| Bulk resins | | |
| WorkBeads 40 Ni-NTA | 25 ml | 40 651 001 |
| WorkBeads 40 Co-NTA | 25 ml | 40 651 401 |
| WorkBeads 40 Cu-NTA | 25 ml | 40 651 301 |
| WorkBeads 40 Zn-NTA | 25 ml | 40 651 504 |
| WorkBeads 40 Ni-IDA | 25 ml | 40 650 001 |
| WorkBeads 40 Co-IDA | 25 ml | 40 650 401 |
| WorkBeads 40 Cu-IDA | 25 ml | 40 650 301 |
| WorkBeads 40 Zn-IDA | 25 ml | 40 650 501 |
| Accessories | | |
| Column plug male 1/16" | 10 | 70 100 010 |
| Column cap female 1/16" | 10 | 70 100 020 |

1. Other pack sizes can be found in the complete product list on www.bio-works.com

Ordering information

| Product name | Pack size | Article number |
|---|-----------|----------------|
| BabyBio NTA His-tag Screening kit 1 ml ¹ | 1 ml x 4 | 45 700 101 |
| BabyBio NTA His-tag Screening kit 5 ml ¹ | 5 ml x 4 | 45 700 102 |
| BabyBio IDA His-tag Screening kit 1 ml ¹ | 1 ml x 4 | 45 700 001 |
| BabyBio IDA His-tag Screening kit 5 ml ¹ | 5 ml x 4 | 45 700 002 |

1. Includes one column each charged with Ni²⁺, Co²⁺, Cu²⁺ or Zn²⁺

Orders: sales@bio-works.com or contact your local distributor.

For more information about local distributor and products please visit www.bio-works.com or contact us at info@bio-works.com



Bio-Works
Virdings allé 18
754 50 Uppsala
Sweden